

# **MANUAL: Materials and Mechanics Practice (PHY109P)**

## **B. Tech: July---November 2017**

**IIITD&M Kancheepuram**  
**Chennai – 600 127**



# LIST OF EXPERIMENTS:

## **CYCLE – I**

1. TORSIONAL PENDULUM
2. BAR PENDULUM
3. CREEP TEST
4. TENSILE TEST
5. FRICTION

## **CYCLE – II**

6. STRAIN GAUGE
7. THREE POINT BEND TEST
8. BUCKLING TEST
9. MICROSTRUCTURE
10. HARDNESS TEST AND  
TORQUE MEASUREMNET

## EX. NO.1: TORSIONAL PENDULUM

### OBJECTIVE:

To find the modulus of rigidity ( $\eta$ ) and torsional rigidity ( $C$ ) of the given string.

### APPARATUS REQUIRED:

Circular or rectangular discs suspended from a point using a metal wire about an axis passing through the middle of plate having largest area, two similar weights, a stop watch and meter scale.

### BRIEF DISCUSSION AND RELEVANT FORMULA:

For an object under torsion torsional, restoring torque is proportional to the angular displacements. In practice, this will be true only for **small** angular displacements ' $\theta$ '.

$$\tau = I_{tot} d^2\theta/dt^2 = -C\theta \quad (1)$$

Again from Fig.2, we get the rigidity modulus  $\eta = \frac{\left(\frac{dF}{\pi r^2}\right)}{\left(\frac{r\theta}{L}\right)} = \frac{l(rdF)}{\pi r^4\theta} = \frac{l\tau}{\pi r^4\theta}$  (2)

We can write the magnitude of  $C$  from (1) and (2),  $C = \frac{\pi \eta r^4}{2l}$  (3)

The equation gives the time period of torsional oscillations of the system as,

$$T = 2\pi\sqrt{\frac{I_{tot}}{C}} = 2\pi\sqrt{\frac{I_o + 2I_s + 2m_s x^2}{C}} \quad (4)$$

$I_o$  = moment of inertia of large disc (without any mass added to it)

$I_s$  = moment of inertia of weight - about an axis passing through its/their own center of gravity, parallel to its/their length.

$m_s$  = mass of each solid cylinder (weight)

$x$  = distance of each weight from axis of suspension.

$C$  = torsional rigidity of suspension wire - (couple per unit twist).

Fig.1





Fig.2



Squaring (4), we get

$$T^2 = \frac{4\pi^2}{C} [I_o + 2I_s + 2m_s x^2] = \frac{8\pi^2 m_s x^2}{C} + \frac{4\pi^2 (I_o + 2I_s)}{C} \quad (5)$$


A graph is plotted between  $T^2$  and  $x^2$ . Fit a best - fit straight line to the data.

$$\text{For the straight line, } slope = \frac{\Delta T^2}{\Delta x^2} = \frac{8\pi^2 m_s}{C} \text{ or } C = 8\pi^2 m_s \frac{\Delta x^2}{\Delta T^2} \quad (6)$$

Knowing, C from (6) one can calculate the rigidity modulus ( $\eta$ ) from (3).

**Note that:**

- The suspension wire should be free from kinks.
- The system should be horizontal always.
- Solid weights must be identical.
- Oscillations should be purely rotational.
- The wire should not be twisted beyond elastic limits.
- The time period should be noted carefully taking the average of about ten periods.
- Make sure that the angular displacements are small.



The student must then compute,

Couple per unit twist of suspension wire C

**Determination of  $I_s$  and  $I_o$ :**  $I_s$  and  $I_o$  can be obtained by measuring corresponding time periods as mentioned below.

Let  $T_o$  = time period of oscillation of large disc alone.

$T_x$  = time period of oscillations of large disc with masses at distance x.

$T_{x0}$  = time period of oscillation of large disc with masses at  $x = 0$ , (obtained from graph – y intercept).

Putting these time periods in (5), we will get a set of equations which will give us the following expressions:

$$I_0 = 2m_s x^2 \frac{T_0^2}{T_x^2 - T_{x0}^2} = \frac{2m_s T_0^2}{\text{slope}}$$

$$I_s = m_s x^2 \frac{T_{x0}^2 - T_0^2}{T_x^2 - T_{x0}^2} = \frac{m_s (T_{x0}^2 - T_0^2)}{\text{slope}}$$

#### OBSERVATIONS:

Mass of identical weights added,  $m_s = \dots\dots\dots$  kg

Time period of disc with no masses  $T_0 = \dots\dots\dots$  sec (averaged over 10 cycles)

Y-intercept from graph  $T_{x0} = \dots\dots\dots$  sec

**TABLE 1: MEASURE THE TIME PERIOD**

Distance of weight from axis of twist $x$ (cm)	$x^2$ (cm <sup>2</sup> )	Number of oscillations, $n$	Time taken $t$ (sec)	Period $T = t / n$ (sec)	$T^2$ (sec <sup>2</sup> )
$x_1$					
$x_2$					
$x_3$					
$x_4$					

**RESULT:** 1. Rigidity modulus of the wire,  $\eta =$

2. Couple per unit twist,  $C =$

3. Moment of Inertia of the disc,  $I_0 =$

## EX.NO. 2: BAR PENDULUM

### OBJECTIVE:

1. To determine the acceleration due to gravity ( $g$ ) using a bar pendulum.
2. To verify that there are two pivot points on either side of the centre of gravity (C.G.) about which the time period is the same.
3. To determine the radius of gyration of a bar pendulum by plotting a graph of time period of oscillation against the distance of the point of suspension from C.G.
4. To determine the length of the equivalent simple pendulum.

### APPARATUS REQUIRED:

Bar Pendulum, Small metal wedge, Spirit level, Stop watch, Meter rod


### THEORY:

A bar pendulum is the simplest form of compound pendulum. It is in the form of a rectangular bar (with its length much larger than the breadth and the thickness) with holes (for fixing the knife edges) drilled along its length at equal separation.

If a bar pendulum of mass  $M$  oscillates with a very small amplitude  $\vartheta$  about a horizontal axis passing through it, then its angular acceleration ( $d^2\vartheta/dt^2$ ) is proportional to the angular displacement  $\vartheta$ . The motion is **simple harmonic** and the time period  $T$  is given by

$$T = 2\pi \sqrt{\frac{I}{Mgl}} \quad (1)$$

where  $I$  denotes the **moment of inertia** of the pendulum about the horizontal axis through its **center of suspension** and  $l$  is the distance between the center of suspension and C.G. of the



Photograph of a typical bar pendulum

Pendulum. According to the theorem of parallel axes, if  $I_G$  is the moment of inertia of the pendulum about an axis through C.G., then the moment of inertia  $I$  about a parallel axis at a distance  $l$  from C.G. is given by

$$\begin{aligned} I &= I_G + Ml^2 \\ &= Mk^2 + Ml^2 \end{aligned} \quad (2)$$

Where,  $k$  is the **radius of gyration** of the pendulum about the axis through C.G. Using Equation (2) in Equation (1), we get

$$\begin{aligned} T &= 2\pi \sqrt{\frac{Mk^2 + Ml^2}{Mgl}} \\ T &= 2\pi \sqrt{\frac{k^2 + l^2}{gl}} = 2\pi \sqrt{\frac{k^2/l + l}{g}} = 2\pi \sqrt{\frac{L}{g}} \end{aligned} \quad (3)$$

Where,  $L$  is the length of the equivalent simple pendulum, given by

$$L = \left( \frac{k^2}{l} + l \right) \quad (4)$$

Therefore, from (1) and (4),  $g = 4\pi^2 \frac{L}{T^2}$  (5)

Equation (4) is a quadratic equation for  $l$ , which must have two roots  $l_1$  and  $l_2$  (say) and follows,  $l_1 + l_2 = L$  and  $l_1 l_2 = k^2$  (6)

for any particular value of  $l(l_1)$ , there is a second point on the same side of C.G. and at a distance  $l(k^2/l_1)$  from it, about which the pendulum will have the same time period. And the graph will look like:




Figure 1: Expected variation of time period with distance of the point of suspension from C.G.

### Ferguson's method for determination of g

Using Equations (4) and (5) we get


$$l^2 = \frac{g}{4\pi^2} lT^2 - k^2$$

A graph between  $l^2$  and  $lT^2$  should therefore be a straight line with slope,  $\frac{g}{4\pi^2}$  as shown in (Figure2).

The intercept on the y-axis is  $\rightarrow k^2$ .

Acceleration due to gravity,  $g = 4\pi^2 \times \text{slope}$

Radius of gyration,  $k = \sqrt{(\text{intercept})}$



**Figure2: Expected form of the graph between  $l^2$  and  $lT^2$ .**

#### PROCEDURE:

1. Balance the bar on a sharp wedge and mark the position of its C.G.
2. Fix the knife edges in the outermost holes at either end of the bar pendulum. The knife edges should be horizontal and lie symmetrically with respect to centre of gravity of the bar.
3. Check with spirit level that the glass plates fixed on the suspension wall bracket are horizontal. The support should be rigid.
4. Suspend the pendulum vertically by resting the knife edge at end A of the bar on the glass plate.
5. Displace the bar slightly to one side of the equilibrium position and let it oscillate with the amplitude not exceeding 5 degrees. Make sure that there is no air current in the vicinity of the pendulum.
6. Use the stop watch to measure the time for 20 oscillations. The time should be measured after the pendulum has had a few oscillations and the oscillations have become regular.
7. Measure the distance  $l$  from C.G. to the knife edge.



8. Record the results in Table 1. Repeat the measurement of the time for 20 oscillations and take the mean.
9. Suspend the pendulum on the knife edge of side B and repeat the measurements in steps 5-8 above.
10. Fix the knife edges successively in various holes on each side of C.G. and in each case, measure the time for 30 oscillations and the distance of the knife edges from C.G.

### OBSERVATIONS:

**TABLE 1: MEASUREMENT OF  $T$  AND  $l$**

Least count of stop-watch = .....sec.

S.No.	Side A up					Side B up				
	Time for 20 Oscillations (t)		$t$ (mean)	$T = t/20$ (sec)	$l$ (cm)	Time for 20 Oscillations (t)		$t$ (mean)	$T = t/20$ (sec)	$l$ (cm)
	1	2				1	2			
1										
2										
3										
4										
5										
6										
7										
8										
9										

### Calculations

Plot a graph showing how the time period  $T$  depends on the distance from the center of suspension to C.G. ( $l$ ). **Figure.1** shows the expected variation of time period with distance of the point of suspension from C.G.

### Acceleration due to gravity ( $g$ )

Draw horizontal lines on the graph corresponding to a period,  $T_1$  as shown in (Figure 1). For the line ABCDE

**Radius of gyration ( $k$ ) and  $g$  calculation:**

$$\text{Let } l_1 = \frac{1}{2}(AC + CE) = \frac{1}{2}AE,$$

$$\text{and } l_2 = \frac{1}{2}(BC + CD) = \frac{1}{2}BD. \quad (7)$$

$$\text{And from (6), } L = \frac{1}{2}(AE + BD)$$

Hence, using this  $T_1$  and  $L$  in the formula for  $g$  (5) we get,

$$g = \text{.....cm/sec}^2.$$

The radius of gyration (6) can be evaluated using using the expression

$$k = \sqrt{l_1 l_2} = \dots\dots\dots \text{cm.}$$

Repeat it for another line A'B'C'D'E' (say for period  $T_2$ ) and calculate the mean values for g and k.

If M is the mass of the bar pendulum, the moment of inertia of the bar pendulum is obtained using the equation

$$I = Mk^2$$

For **Ferguson's method**, fill up the following table to evaluate  $I^2$  and  $IT^2$  corresponding to all the measurements recorded in Table 1.

**TABLE 2: CALCULATED VALUES OF ( $I^2$ ) AND ( $IT^2$ ):**

S.No.	Side A up		Side B up		Mean values	
	$I^2$ (cm <sup>2</sup> )	$IT^2$ (cm sec <sup>2</sup> )	$I^2$ (cm <sup>2</sup> )	$IT^2$ (cm sec <sup>2</sup> )	$I^2$ (cm <sup>2</sup> )	$IT^2$ (cm sec <sup>2</sup> )
1						
2						
3						
4						
5						
6						
7						
8						
9						

Plot a graph of  $I^2$  against  $IT^2$  (as shown in Figure.2) and determine the values of the slope and the intercept on the  $I^2$  axis.

#### Interference Obtained from Graph:

Slope of the graph = ..... cm/sec<sup>2</sup>.

Intercept = ..... cm<sup>2</sup>.

Acceleration due to gravity  $g = 4\pi^2 \times \text{slope} = \dots\dots\dots \text{cm/sec}^2$ .

Radius of gyration,  $k = \sqrt{(\text{intercept})} \dots\dots \text{cm}^2$ .

### EX.NO. 3: CREEP TEST

AIM: To study the behavior of the material under long-term constant loading and to determine the minimum creep rate (in Stage II).

APPARATUS REQUIRED: Creep Test-Rig, Micrometer, Vernier caliper, Material for testing, Weights

FORMULA: A creep test involves a tensile specimen under a constant load maintained at a constant temperature. Measurements of strain are then recorded over a period of time.

$$\text{Creep strain} = \text{Elongation} / L_0$$

$$\text{Applied stress} = \text{Load/Area ( N/mm}^2\text{)},$$

where,  $L_0$  and  $A_0$  are, respectively, the initial length and area of the specimen

$$\text{Creep modulus} = \text{Applied Stress} / \text{Creep Strain}$$

#### PROCEDURE:

1. Make the sample by ASTM standard
2. Measure the length and thickness of the given test specimen.
3. Set the specimen in between the jaws.
4. Move the dial gauge and make contact with a plate attached to the movable jaw and the make the dial gauge reading as zero.
5. Hold the plunger and place the weight gradually and then release the plunger.
6. Note down the dial gauge reading for every 30 seconds till the given specimen breaks.
7. Tabulate the readings and determine the creep strain and creep modulus for each dial gauge reading.




TABLE 1: TO FIND THE CREEP MODULUS


Material Response for Stress (I) at 3.0 Kg:

S.NO	TIME (s)	ELONGATION (mm)	CREEP STRAIN	CREEP MODULUS (N/mm <sup>2</sup> )

- Repeat the same experiment with a different load for Stress(II) at 3.5kg.

GRAPH:

With common X & Y axis draw the creep curves (time Vs strain) for the material at different loads



RESULTS:

Initial Creep strain for the material with stress  $\sigma_1 =$

Minimum creep rate for the material with stress  $\sigma_1 =$   $s^{-1}$

INFERENCE:

Compare the results for different initial load and draw the schematic representation of the effect of stress on creep curves at a constant temperature.

#### EX. 4: TENSILE TEST

AIM: To study the response of the given specimens subjected to tensile load.

DESCRIPTION: The engineering tension test is widely used to provide basic design information on the strength of the materials and as an accepted test for the specification of the materials. In the tension test, a specimen is subjected to a continually increasing uniaxial tensile force while simultaneous observations are made of the elongation of the specimen.

APPARATUS/ INSTRUMENT REQUIRED: INSTRON tensile testing machine, Capacity-2 kN, Vernier caliper and scale, and Test specimens- As per ASTM standards

#### PROCEDURE:

1. Measure and record the initial dimension of the specimen (gauge length- $L_0$ , width  $w_0$ , thickness  $t_0$ , cross section area  $A_0 = w_0 \times t_0$ ).
2. Fix the test specimen between fixed and movable jaws of the machine.
3. Reset the load to zero.
4. Operate the machine till the specimen fractures.
5. Measure and record the final configuration of the specimen (gauge length  $L_f$ , width  $w_f$ , thickness  $t_f$ , cross section area  $A_f = w_f \times t_f$ ).
6. Repeat the experiment for different strain rate (rate of loading).
7. Using the data acquired by the system, construct the stress-strain curves and find the various parameters as listed in the calculation.

#### OBSERVATION:

Sl. No	Material, Strain rate & Load	Dimensions (mm)	Fracture dimension (mm)
1	<u>Aluminium</u>  Strain rate:  Load:	$L_0 =$  $t_0 =$  $w_0 =$  $A_0 =$	$L_f =$  $t_f =$  $w_f =$  $A_f =$
2	<u>Nylon</u>  Strain rate:  Load:	$L_0 =$  $t_0 =$  $w_0 =$  $A_0 =$	$L_f =$  $t_f =$  $w_f =$  $A_f =$

#### CALCULATIONS:

1.  $A_0 = \quad mm^2$

2.  $A_f = \quad mm^2$

3. Ultimate tensile strength,  $S_u = \frac{P_{max}}{A_0} = \quad N / mm^2$

Where,  $P_{max}$  is the maximum load

4. Yield strength,  $S_0 = \quad N / mm^2$  (obtain from graph)

Note: Yield strength is the stress required to produce a small specified amount of plastic deformation. The usual definition of this property is the offset yield strength determined by the stress corresponding to the intersection of the stress-strain curve and a line parallel to the elastic part of the curve offset by a strain of 0.2%.

5. Breaking stress,  $S_f = \frac{P_f}{A_0} = \quad N / mm^2$

where,  $P_f$  is the breaking/fracture load (load at the occurrence of fracture)

6. Strain,  $e_f = \frac{L_f - L_0}{L_0} =$

7. Reduction in area at fracture,  $q = \frac{A_0 - A_f}{A_0} =$


8. Modulus of elasticity,  $E = \text{Slope of initial linear portion of the curve,} \quad N / mm^2$

9. Resilience,  $U_R = \frac{S_0^2}{2E} = \quad N / mm^2$

10. Toughness,  $U_T = \frac{S_0 + S_u}{2} e_f = \quad N / mm^2$

#### PLOTS:

1. Engineering stress Vs Engineering strain
2. True stress Vs True strain



#### INFERENCE:

1. Compare the results and state which material has high strength, toughness, ductility, and stiffness.
2. State the effect of strain rate in material response.



### EX.5: DETERMINATION OF COEFFICIENT OF STATIC FRICTION

**OBJECTIVE:** To measure the static coefficient of friction for several combinations of material surfaces.

**APPARATUS REQUIRED:** Inclined plane, Metal block, pull-push meter, set of weights and materials with different surfaces.

**FORMULA:** Coefficient of static friction,  $\mu_s = \frac{f_s}{N}$

Angle of static friction,  $\phi_s = \tan^{-1}(\mu_s)$


where,

$f_s$  – Maximum static friction force (N)

$N$  – Normal force applied (N)



**THEORY:**

The coefficient of static friction ( $\mu_s$ ) can be measured experimentally for an object placed on a flat surface and pulled by a known force. The friction force depends on the coefficient of static friction and the Normal force ( $N$ ) on the object from the surface. When the object just begins to slide, the friction force will attain its maximum value as shown in Figure 1a.



**Figure 1a**

The forces acting on a body kept on an inclined plane are shown in Figure 1.  $W$  is the weight of the body ( $W = mg$ ),  $N$  is the normal force from the plane and  $f$  is the frictional force. Generally, this situation is analyzed by resolving the forces into components parallel and perpendicular to the plane, as shown in Figure 1.



**Figure 1.**

The weight  $W$  is resolved into two components acting along the plane  $W_P$  and normal to plane  $W_N$  which is balanced by frictional force  $f$  and the normal force  $N$ , respectively while the body is stationary. When motion is impending, the friction force  $f$  attains its maximum value  $f_s$ .  
 $f_s = W_P = W \sin \theta_s$

The angle at which the motion is impending is called the 'Angle of repose'.

$$\mu_s = \frac{f_s}{N} = \frac{W \sin \theta_s}{W \cos \theta_s} \rightarrow \theta_s = \tan^{-1}(\mu_s),$$

The angle of repose is equal to the angle of static friction.

#### Procedure

1. Place the metal block on the rubber sheet given.
2. Place some weight at the centre of the metal block and pull it horizontally using the pull-push meter.
3. Note the reading shown on the pull-push meter as motion is impending.
4. Keep the metal block on an inclined wooden plane, whose initial inclination does not exceed  $10^\circ$ , with the rubber sheet between metal block and wooden plane.
5. Slowly increase the angle of the wooden plane.
6. Note the inclination at which motion of metal block is impending, i.e. the angle of repose for the given condition.
7. Increase the load on the metal block and repeat the procedure from step 1.
8. Above experiment can be repeated for different material surfaces.

TABLE 1: TO FIND THE COEFFICIENT OF STATIC FRICTION ( $\mu_s$ ) ON HORIZONTAL PLANE

Trial No.	Surface Type	Total Weight of metal block, $W = N$	Max. friction force, i.e. Pull-Push meter reading, $f_s$	Co-efficient of static friction, $\mu_s = f_s / N$	Angle of static friction, $\phi_s = \tan^{-1}(\mu_s)$
1	-----				
2					
3					

4					
5					

\*\* Repeat this for other given surfaces also

TABLE 2: TO FIND THE COEFFICIENT OF STATIC FRICTION ( $\mu_s$ ) ON INCLINED PLANE

Trial No.	Surface Type	Length ( $l$ )	Height ( $h$ )	Angle of repose, ( $\theta_s$ )	$\mu_s = \tan(\theta_s)$
1	-----				
2					
3					

\*\* Repeat this for other given surfaces

Graph

Plot a graph of the normal force ( $N$ ) and the frictional force ( $f_s$ ) obtained while the metal block was kept on the horizontal plane.

TABLE 3: Comparison of (i) values of  $\mu_s$ , and (ii) the values of  $\theta_s$  and  $\phi_s$ .

Trail No.	Surface Type	Angle of Static Friction ( $\phi_s$ )	Angle of repose, ( $\theta_s$ )	$\mu_s$ (From Horizontal Plane)	$\mu_s$ , (From Inclined Plane)

## EX. 6: YOUNG'S MODULUS OF WOOD USING A STRAIN GAUGE

**AIM:** To determine Young's modulus of a half meter wooden scale using a Strain Gauge.

**APPARATUS REQUIRED:** A half meter scale with two identical strain gauges attached to one end of the scale, one strain gauge at the top and the other at the bottom; another end of the scale is attached to the table with a clamp; a circuit board with appropriate terminals to constitute a Wheatstone bridge network.

### STRAIN GAUGE:

Young's modulus (Y) of the bar (scale) is defined by the ratio of stress (F/A) and tensile strain ( $\Delta L/L$ ),

$$\frac{F/A}{\Delta L/L} = Y \dots \dots \dots (1)$$

where, F is the force applied (Newton), A is the cross-sectional area ( $m^2$ ),  $\Delta L$  is the change in length (m), L is the original change in length (m).

A strain gauge is a transducer whose electrical resistance varies in proportional to the amount of strain in the device. The most widely used gauge is metallic strain gauge which consists of a very fine wire or, more commonly, metallic foil arranged in a grid pattern. The grid pattern maximizes the amount of metallic wire or foil subject to strain in the parallel direction (Fig.1). The cross-sectional area of the grid is minimized to reduce the effect of shear strain and Poisson strain. The grid is bonded to a thin backing, called the carrier, which is attached directly to the test specimen. Therefore, the strain experienced by the test specimen is transferred directly to the strain gauge, which response with a linear change in electrical resistances.




Figure.1

A fundamental parameter of the strain gauge is its sensitivity to strain, expressed quantitatively as the gauge factor ( $\lambda$ ). Gauge factor is defined as the ratio of fractional change in electrical resistance to the fractional change in length (strain).

$$\frac{\Delta R/R}{\Delta L/L} = \lambda \dots \dots \dots (2)$$

The gauge factor ( $\lambda$ ) for metallic strain gauge is typically around 2.

### WHEATSTONE BRIDGE:

Measuring the strain with a strain gauge requires accurate measurement of very small change in resistance and such small changes in  $R$  can be measured with a Wheatstone bridge. A general Wheatstone bridge consists of four resistive arms with an excitation voltage,  $V_{EX}$ , that is applied across the bridge (Figure.2)




Figure 2. Wheatstone Bridge

The output voltage of the bridge,  $V_O$ , will be equal to:

$$V_O = \left[ \frac{R_3}{R_3 + R_4} - \frac{R_2}{R_1 + R_2} \right] \cdot V_{EX} \dots \dots \dots (3)$$

From this equation, it is apparent that when  $R_1/R_2 = R_3/R_4$ , the output voltage  $V_O$  will be zero. Under this condition, the bridge is said to be balanced. Any change in resistance in any arm of the bridge will result in a non-zero output voltage. Therefore, if we replace  $R_4$  in Figure 2 with an active strain gauge, any change in the strain gauge resistance will unbalance the bridge and produce a nonzero output voltage. If the nominal resistance of the strain gauge is designed as  $R_G$ , then the strain induced change in resistance  $\Delta R$ , can be expressed as

$$\Delta R = R_G \cdot \lambda \cdot \text{strain} \dots \dots \dots (4)$$

Assuming that  $R_1 = R_2$  and  $R_3 = R_G$ , the bridge equation above can be rewritten to express  $V_O/V_{EX}$  as a function of strain.

Ideally, we would like the resistance of the strain gauge to change only in response to applied strain. However, strain gauge material, as well as the specimen material on which the gauge is mounted, will also respond to changes in temperature. Strain gauge manufacturers attempt to minimize sensitivity to temperature by processing the gauge material to compensate for the thermal expansion of the specimen material for which the gauge is intended. While compensated gauges reduce the thermal sensitivity, they do not totally remove it. By using two strain gauges in the bridge, the effect of temperature can be further

minimized. For example, in a strain gauge configuration where one gauge is active ( $R_G + \Delta R$ ), and a second gauge is placed transverse to the applied strain. Therefore, the strain has little effect on the second gauge, called the dummy gauge. However, any changes in temperature will affect both gauges in the same way. Because the temperature changes are identical in the two gauges, the ratio of their resistance does not change, the voltage  $V_o$  does not change, and the effects of the temperature change are minimized.

The sensitivity of the bridge to strain can be doubled by making both gauges active in a half bridge configuration. Figure.3 illustrates a bending beam application with one bridge mounted in tension( $R_G + \Delta R$ ) and the other mounted in compression ( $R_G - \Delta R$ ). This half bridge configuration, whose circuit diagram is also illustrated in Figure.3 yields an output voltage that is linear and approximately doubles the output of the quarter bridge.




Figure.3

And in this experiment, we aim to determine Young's modulus of a half meter wooden bar by loading it with a mass of " $m$ " gram. For a beam of rectangular cross section with breadth  $b$  and thickness  $d$ , the moment of inertia  $I$ , is

$$I = b d^3 / 12 \dots \dots \dots (5)$$

The moment of force/restoring couple is  $Y.I / r_c$  where  $r_c$  is the radius of curvature of the bending beam. The Young's modulus is calculated by assuming that at equilibrium, the bending moment is equal to the restoring moment.

#### PROCEDURE:

1. Clamp the beam to the table in such a way that the strain gauges are close to the clamped end. Load and unload the free end of the beam a number of times.
2. Make the connections as given in the circuit diagram (Figure.4)





Figure.4

$P=100\Omega$  resistor,  $S_1=10$  mA current source, DPM= a voltmeter with digital panel.  $R=S$ = strain gauge resistance  $\sim 120\Omega$  with a gauge factor  $\lambda=2.2$ ,  $Q=82\Omega$  (plus the resistance of the rheostats) all in series.

3. Switch ON the constant current source ( $S_1$ ) and DPM.
4. Balance the bridge using the rheostats. At this stage the DPM will read or very nearly zero. Note that, this is done at no load (only with the dead load).
5. Load the beam with a hanger of mass ' $m$ ' gm suspending it as close to the free end of the scale as possible. Note the DPM reading. (Note that as you are about to take a reading the last digit will be changing about the actual steady value. Take at least 10 readings continuously and take the average these ten).
6. Increase the load in steps of  **$m$  (50) gm**, up to  **$5m$**  gm and take the readings each time.
7. Unload the beam from  **$5m$**  down to zero in steps of  **$m$**  gm at a time and note the DPM reading each time.
8. To check reproducibility, repeat all the above process taking readings while loading and unloading in steps of  **$m$**  gm.
9. Draw a graph between  **$m$**  along X axis and unbalanced voltage  **$V$**  along Y axis. Determine the slope of the graph ( **$dV/dm$** )



10. Note the distance between the center of the strain gauges and the point of application of the load (**L**).
11. Measure the breadth of the beam using slide calipers (**b**).
12. Measure the thickness of the beam using a screw gauge (**d**).
13. Young's modulus of the material of the beam, which is nothing but stress to strain ratio, is given by the following expression (Working formula).

$$Y = \frac{6gL\lambda RI}{bd^2[1 + (R/P)] \frac{dV}{dm}} \dots \dots \dots (5)$$

where

- $g$  is the acceleration due to gravity,
- $\lambda$  is the gauge factor (for metal strain gauge  $\lambda = 2.2$ ).
- $I$  is the output current from the source  $S_1$ .
- $R$  is the resistance of strain gauge.
- $\frac{dV}{dm}$  is slope of the  $m$  Vs  $V$  curve

**TABLE 1: Change of voltage with load**

<div style="display: flex; align-items: center;"> <div style="text-align: right; padding-right: 10px;">Load (gram)</div> <div style="text-align: left; padding-left: 10px;">DPM reading (mV)</div> </div>	0m	1m	2m	3m	4m
1) Loading $V_1$					
2) Unloading $V_2$					
3) Mean of $V_1 + V_2$					

**RESULT:**


Thus, Young's modulus of the given wooden scale is,  $Y = \text{-----} N/m^2$ .



### EX.7: THREE POINT BEND TEST

AIM: To determine the modulus of elasticity of metallic and non-metallic materials in bending and to verify Maxwell's reciprocal theorem.

APPARATUS REQUIRED: Three-point bend test apparatus, Micrometer, Vernier caliper, Weights



#### FORMULA:

$$E = \frac{WL^3}{48Iy_{\max}}$$

Where,

$E$  – Modulus of elasticity in bending in  $N / mm^2$ .

$W$  – Load in  $N$ .

$L$  – Span length in  $mm$ .

$I$  – Moment of inertia of cross-section in  $mm^4$ .

For Rectangular bar,  $I = \frac{bd^3}{12}$

Where,  $b$  – breadth,  $d$  – thickness of the specimen in  $mm$

For Cylindrical rod,  $I = \frac{\pi d^4}{64}$

Where,  $d$  – diameter of the specimen in  $mm$

and

$y_{\max}$  - Maximum deflection under point load in  $mm$ .

#### PROCEDURE:

1. Measure the dimensions of the test specimen.
2. Fix the support at  $L$  distance apart.
3. Place the specimen on the support.
4. Rotate the dial gauge outer frame and ensure zero setting.
5. Place the loading table exactly at the mid-span of the specimen.

6. Place the weight on the loading pan and note down the total load and deflection from the dial gauge.
7. Increase the load in steps and measure the deflection.
8. Vary the span length and repeat the experiment for each specimen.
9. Calculate the modulus of elasticity in bending.
10. Repeat the experiment for different specimens.

OBSERVATION:

**TABLE 1: TO FIND THE BREADTH OF THE SPECIMEN: LC= \_\_\_\_\_**

S. No	MSR (mm)	VSC (div)	VSR = (VSC*LC) (mm)	TR = MSR+VSR (mm)
Mean				

**TABLE 2: TO FIND THE THICKNESS OF THE SPECIMEN : LC= \_\_\_\_\_**

S. No	MSR (mm)	VSC (div)	VSR = (VSC*LC) (mm)	TR = MSR+VSR (mm)
Mean				

**TABLE 3: TO FIND THE DIAMETER OF THE SPECIMEN USING SCREW GAUGE: LC= \_\_\_\_\_**

S. No	PSR (mm)	HSC (div)	HSR = (VSC*LC) (mm)	TR = PSR+HSR (mm)	CR= TR±ZE (mm)
Mean					

**TABLE 4: TO FIND THE MODULUS OF ELASTICITY FOR DIFFERENT MATERIALS**


Span Length (mm)	Load (N)	Deflection (mm)	Modulus of elasticity in bending (N/mm <sup>2</sup> )	Avg. Modulus of elasticity in bending (N/mm <sup>2</sup> )
Specimen-1: _____				
Moment of Inertia, I=				
Specimen-2:				
Moment of Inertia, I=				

Specimen-3:				
Moment of Inertia, I=				

**CALCULATIONS:** Use formulae ---


$$I = \frac{bd^3}{12} \quad I = \frac{\pi d^4}{64} \quad E = \frac{WL^3}{48I_{y_{\max}}}$$

**GRAPH:** Draw the graph of load(P) vs deflection(y) for different specimens using common X & Y axis.



### MAXWELL'S RECIPROCAL THEOREM:

The Maxwell's law of reciprocal deflection states that, if  $a$  and  $b$  are two points in a loaded beam, then the deflection ( $y$ ) at  $b$  due to application of unit load  $P$  at  $a$  is exactly equal to the deflection ( $y$ ) at  $a$  due to the same unit load  $P$  applied at  $b$ . (Similarly for the points  $c$  and  $d$  also)



### VERIFICATION OF MAXWELL'S RECIPROCAL THEOREM:

1. Measure the dimensions of the test specimen.
2. Fix the support at  $L$  distance apart.
3. Divide the span length into six equal segments and place the load ( $P$ ) at a distance ' $a$ ' from left support.
4. Measure & tabulate the deflection at beam at a distance ' $a$ ', ' $b$ ', ' $c$ ' and, ' $d$ ' from left support.
5. Conduct the experiment for different loads & tabulate the deflection.

**TABLE 5: VERIFICATION MAXWELL'S RECIPROCAL THEOREM**

<b>Specimen:</b>		<b>Dimension:</b> $l = \text{-----}; b = \text{-----}; t = \text{-----}$		<b>Span length:</b>	
<b>Load location (from left support)</b>	<b>Load</b>	<b>Deflection</b>			
		<i>a</i>	<i>c</i>	<i>b</i>	<i>d</i>
<i>a</i>					
<i>c</i>					
<i>b</i>					
<i>d</i>					

Repeat this TABLE for different materials with different dimensions.

INFERENCES:


1. State the reason for variation in  $E_{avg}$  value in different span length
2. State how the cross section of the beam affects the deflection

## EX.08: BUCKLING OF COLUMNS EXPERIMENT

**OBJECTIVE:** To observe the buckling behavior of columns and estimate their buckling loads for different end conditions

**APPARATUS REQUIRED:** Structural testing frame set-up, columns made of different materials, weights, traveling microscope, vernier calipers and screw gauge.

**THEORY:** Buckling is a stability problem. If a column (rod) is subjected to longitudinal forces, as shown in Figure below, it can fail in two ways: 1) it can be plasticized and flattened if its admissible compressive strain is exceeded called compressive failure or 2) it is possible that it will suddenly shift to one side and buckle before attaining the admissible compressive strain. The second way of failure is known as buckling. The shape of the rod is the factor that determines which of the two cases of failure will occur. A slender, thin rod is more likely to buckle than a thick, stout rod.



The critical load  $P_{CR}$  above which buckling can occur is dependent on length, diameter (width), end conditions of the rod and the material used. According to Euler's theorem of buckling,

$$\text{Critical Load, } P_{CR} = \pi^2 E I / (L_e)^2 \quad (1) \text{ where,}$$

$E$  – Elasticity Modulus of material ( $\text{N/mm}^2$ )

For Polypropylene, PP material =  $1.5 \cdot 10^9 \text{ N/m}^2$

For Stainless Steel =  $180 \cdot 10^9 \text{ N/m}^2$

$I$  – Least Moment of Inertia of cross section of column ( $\text{mm}^4$ )

For Rectangular bar,  $I = \frac{bd^3}{12}$ , with,

$b$  – Breadth of the specimen (mm)


$d$  – Thickness of the specimen (mm)

$L$  – Actual Length of column (mm)

$L_e$  – Effective Length:


$L_e = 2L$  for Fixed-Free;

$= 1.0 L$  for both Pinned  
 $= 0.7 L$  for Fixed-Pinned  
 $= 0.5 L$  for both Fixed.



Euler derived the equation (1) for columns with no consideration for lateral forces. However, if lateral forces are taken into consideration the value of critical load remains approximately the same. In deriving equation (1), initial lateral deflection of the column has also been taken zero. If the initial deflection of the mid point of the column is 'a' and instantaneous mid point deflection is  $\delta$ , then the critical load can be expressed in terms of its load and deflection of mid-point as (Southwell plot);

$$\delta = \frac{a}{\frac{P_{CR}}{P} - 1} \Rightarrow \frac{\delta}{P} = \left( \frac{1}{P_{CR}} \right) \delta + \left( \frac{a}{P_{CR}} \right) \quad (2)$$



PROCEDURE



1. Measure the length, breadth, and thickness of the column which is to be used in the experiment.
2. Note the Elasticity Modulus of the material given to you.
3. Mount the column providing the required end condition (You may use fixed-free end condition at first).
4. Focus the telescope on the column.
5. Gradually apply a load in small steps on the column using weights on top of the frame set-up and note down the load and the corresponding lateral deflection of the column using dial gauge.
6. Continue step 4 till the lateral deflection is noticeable, i.e. column is no longer straight, and then remove the load.
7. Change the end condition to fixed-fixed, pinned-pinned, pinned-fixed, and repeat the experiment from step 3.
8. Repeat the experiment with columns made of different materials and different cross sections.

**TABLE-1:** Measurement of deflection with a change in load.

**MATERIAL:**

Sl.No	End conditions	Length of column (mm)	Applied Load, P (kgf)	Lateral Deflection, $\delta$ (mm)	Ratio $\delta/P$ , (mm/kgf)
	<b>Fixed-Fixed</b>				
	<b>Fixed-Pinned</b>				

	<b>Pinned-Pinned</b>				

- Repeat this experiment with different material.


**Table-2: TO FIND THE BREADTH OF THE SPECIMEN:** LC= \_\_\_\_\_

S. No	MSR (mm)	VSC (div)	VSR = (VSC*LC) (mm)	TR = MSR+VSR (mm)
<b>Mean</b>				

**Table-3: TO FIND THE THICKNESS OF THE SPECIMEN:** LC= \_\_\_\_\_

S. No	MSR (mm)	VSC (div)	VSR = (VSC*LC) (mm)	TR = MSR+VSR (mm)
<b>Mean</b>				

\*Plot a graph between the measured lateral deflection (abscissa) and the ratio ( $\delta/P$ ). The inverse of the slope of the curve gives the experimental critical load (using Data from Table1).



**TABLE 4: ERROR calculation.**

Sl.No	End conditions	Length of column (mm)	Effective length of column(mm)	Theoretical Critical Load, $P_{CR}$ (N)	Experimental Critical Load $P_{CR}$ (N)	Error (%)

**Caution:**

DO NOT overload the column beyond the critical value (and/or DO NOT attempt for high lateral deflections) as it may result in permanent plastic deformation to the column and/or breaking of columns.

**RESULT/conclusion:** for example remark about theoretical and experimental values of  $P_{CR}$ .

### EX.9: MICROSTRUCTURE PRACTICE

OBJECTIVE: To see and study the microstructure for cast iron and steel 304.

APPARATUS REQUIRED: Sample, Grinding machine, Emery sheets with different grades, Metallurgical Inverted Microscope, Polishing machine, Alumina powder, Etchant, Cotton.

THEORY: Internal structure of a material is viewed on a Microscopic scale. Microstructure refers to the fine surface structure of a pure metal or alloy, as revealed by magnifications of 25X or greater.

The microstructure of a material (such as metals, polymers, ceramics or composites) can strongly influence physical properties such as strength, toughness, ductility, hardness, corrosion resistance, high/low-temperature behavior or wear resistance.

PROCEDURE:

Steps: - Moulding

1. Make the given material into small pieces for testing
2. Prepare the molding to hold the sample.

Moulding Preparation: Mix the cold setting compounds into a colloidal form, use both liquid and powder in a proper mixture.

3. Place the small piece of test material at the bottom of the mould ring, and then pour the prepared colloidal mixture on it.
4. Leave that for 15 – 20 minutes to get harden moulding
5. The sample is ready for the test procedure.

### Steps: - Testing

1. Remove the roughness and micro burs of the sample with the help of belt grinding machine
2. Polish the sample with various grades of emery sheet like 320, 400, 600, 800 and higher grit abrasives, until the desired finish is achieved.
3. Use the alumina powder on disc polishing machine, to get the smooth and impressive surface on the sample.
4. Clean the specimen with cotton
5. Apply the etchant – NITAL (98% Ethanol & 2% Nitric acid)
6. Check the microstructure and draw the structure you see for each sample and comment on it.




Fig: Sample work piece to view microstructure




Fig.1 – Mild Steel

Grains

Grain

Boundary




Fig.2 – Cast Iron

### EX.10a: TORQUE MEASUREMENTS

OBJECTIVE: To understand the working principle of strain gauge and to measure the torque using strain gauge

APPARATUS REQUIRED: Torque measurement kit with digital indicator, Weights, Meter Scale

FORMULA:

$$\tau = r \times F \quad (1)$$

where

$\tau$  – Torque measured (kg-m)

$r$  – radius from the origin point (m)

$F$  – force applied (kg)

PROCEDURE: TORQUE MEASUREMENT


1. Switch on the instrument.
2. Using coarse and fine adjusting knobs set the reading in to zero.
3. Place the weight on the load pan and note down the torque values for a particular length.
4. Note down the readings during successive loading and unloading and take the average value for the particular load.
5. Repeat this for other length values.
6. Compare the strain gauge readings with theoretical torque.

TABLE 1: Variation of torque with load.

S.NO	TORQUE ARM LENGTH (m)	LOAD (g)	STRAIN GAUGE READING (kg-m)			THEORETICAL TORQUE (kg-m)
			LOADING	UNLOADING	AVERAGE	
1	1					
2	0.75					
3	0.5					
4	0.25					

GRAPH:

- Plot the graph of Load Vs Strain gauge measured torque & theoretical torque on common X & Y axis graph.



RESULT:

The working principle of the strain gauge is understood through measuring strain and torque using the strain gauge in a Wheatstone bridge circuit.

INFERENCE:

As for example “State the reason for the variation in theoretical & experimental strain & torque values.”



## EX.10b: HARDNESS TEST

**OBJECTIVE:** To find the Hardness number of different materials

**APPARATUS REQUIRED:** Digital Rockwell Hardness Test Machine, Indenters like Diamond tip and steel ball with various diameters, Different materials for testing

**THEORY:** Hardness refers to various properties of matter in the solid phase that gives it high resistance to various kinds of shape change when force is applied. The hard matter is contrasted with the matter. Macroscopic hardness is generally characterized by strong intermolecular bonds. However, the behavior of solid materials under force is complex, resulting in several different scientific definitions of what might be called "hardness" in everyday usage.

In materials science, there are three principal operational definitions of hardness:

- Scratch hardness: Resistance to fracture or plastic (permanent) deformation due to friction from a sharp object
- Rebound hardness: Height of the bounce of an object dropped on the material, related to elasticity.
- Indentation hardness: Resistance to plastic (permanent) deformation due to a constant load from a sharp object

The equation based definition of hardness is the pressure applied over the projected contact area between the indenter and the material being tested. As a result hardness values are typically reported in units of pressure, although this is only a "true" pressure if the indenter and surface interface is perfectly flat.

The Rockwell scale is a hardness scale based on the indentation hardness of a material. The Rockwell test determines the hardness by measuring the depth of penetration of an indenter under a large load compared to the penetration made by a preload. There are different scales, which are denoted by a single letter (say X), that use different loads or indenters. The result, which is a dimensionless number, is noted by HRX where X is the scale letter.

The Hardness number is,  $= HR X (NUMBER)$ . Where,

HR – Rockwell Hardness, X – Scale corresponds, Number–numerical digit for the particular material

When testing metals, indentation hardness correlates linearly with tensile strength. This important relation permits economically important non-destructive testing of bulk metal

deliveries with lightweight, even portable equipment, such as hand-held Rockwell hardness testers.

### Rockwell Hardness Test



Fig.1 Hardness testing machine




Fig.2 Hardness testing process

The Rockwell hardness test method consists of indenting the test material with a diamond cone or hardened steel ball indenter. The indenter is forced into the test material under a preliminary minor load  $F_0$  (Fig. 2a) usually 10 kg. When equilibrium is reached, an indicating device, which follows the movements of the indenter and so responds to changes in depth of penetration of the indenter, is set to a datum position. While the preliminary minor load is still applied, an additional major load is applied resulting an increase in penetration (Fig. 2b). When equilibrium is reached, the additional major load is removed but the preliminary minor load is still maintained. Removal of the additional major load allows a partial recovery, so reducing the depth of penetration (Fig.2c). The permanent increase in depth of penetration, resulting from the application and removal of the additional major load is used to calculate the Rockwell hardness number.

$$HR = E - e$$

$F_0$  = preliminary minor load in kgf

$F_1$  = additional major load in kgf

$F$  = total load in kgf

$e$  = permanent increase in depth of penetration due to major load  $F_1$  measured in units of

0.002 mm

$E$  = a constant depending on form of indenter: 100 units for diamond indenter, 130 units for steel ball indenter

HR = Rockwell hardness number

$D$  = diameter of steel ball

The Rockwell Hardness Scales —The various Rockwell scales and their applications are shown in the Table 1. The type of penetrator and load used with each are shown in Table 2, which give comparative hardness values for different hardness scales.

**TABLE .1**

Scale	Testing Application
A	For tungsten carbide and other extremely hard materials. Also for thin, hard sheets.
B	For materials of medium hardness such as low- and medium-carbon steels in the annealed condition.
C	For materials harder than Rockwell B-100.
D	Where a somewhat lighter load is desired than on the C scale, as on case-hardened pieces.
E	For very soft materials such as bearing metals.
F	Same as the E scale but using a $\frac{1}{16}$ -inch ball.
G	For metals harder than tested on the B scale.
H & K	For softer metals.
15-N; 30-N; 45-N	Where a shallow impression or a small area is desired. For hardened steel and hard alloys.
15-T; 30-T; 45-T	Where a shallow impression or a small area is desired for materials softer than hardened steel.

**TABLE 2. (VALUES ARE OBTAINED FROM THE MANUFACTURER)**

	Scale Symbol	Indenter Type (Ball dimensions indicate diameter.)	Preliminary Force N (kgf)	Total Force N (kgf)	Typical Applications
Regular Rockwell Scales	A	Spheroconical Diamond	98.07 (10)	588.4 (60)	Cemented carbides, thin steel, and shallow case hardened steel.
	B	Ball - 1.588 mm (1/16 in.)	98.07 (10)	980.7 (100)	Copper alloys, soft steels, aluminum alloys, malleable iron, etc.
	C	Spheroconical Diamond	98.07 (10)	1471 (150)	Steel, hard cast irons, pearlitic malleable iron, titanium, deep case hardened steel, and other materials harder than HRB 100.
	D	Spheroconical Diamond	98.07 (10)	980.7 (100)	Thin steel and medium case hardened steel, and pearlitic malleable iron
	E	Ball - 3.175 mm (1/8 in.)	98.07 (10)	980.7 (100)	Cast iron, aluminum and magnesium alloys, and bearing metals
	F	Ball - 1.588 mm (1/16 in.)	98.07 (10)	588.4 (60)	Annealed copper alloys, and thin soft sheet metals.
	G	Ball - 1.588 mm (1/16 in.)	98.07 (10)	1471 (150)	Malleable irons, copper-nickel-zinc and cupro-nickel alloys.
	H	Ball - 3.175 mm (1/8 in.)	98.07 (10)	588.4 (60)	Aluminum, zinc, and lead.
	K	Ball - 3.175 mm (1/8 in.)	98.07 (10)	1471 (150)	Bearing metals and other very soft or thin materials. Use smallest ball and heaviest load that does not give anvil effect.
	L	Ball - 6.350 mm (1/4 in.)	98.07 (10)	588.4 (60)	
	M	Ball - 6.350 mm (1/4 in.)	98.07 (10)	980.7 (100)	
	P	Ball - 6.350 mm (1/4 in.)	98.07 (10)	1471 (150)	
	R	Ball - 12.70 mm (1/2 in.)	98.07 (10)	588.4 (60)	
	S	Ball - 12.70 mm (1/2 in.)	98.07 (10)	980.7 (100)	
	V	Ball - 12.70 mm (1/2 in.)	98.07 (10)	1471 (150)	

**TABLE 3: TO FIND THE HARDNESS VALUE FOR DIFFERENT MATERIALS**

S.No	Material	Indenter	Scale	Load	Hardness Number

**RESULT:** The Harness numbers for the given materials have been found by using Digital Rockwell Hardness Test machine.

